

## PRODUCT INFORMATION

---

**Product Name :** 6× BPB loading Dye

<b>Code No.</b>	<b>Size</b>
DM210	1 mL × 4
DM212	1 mL × 2

This product is research use only.

**Description :** 6× BPB loading Dye is useful for loading DNA samples into wells on agarose gel electrophoresis. It contains 36 mM EDTA • 2Na, 40 % glycerol, 0.075 % bromphenol blue and 10 mM Tris-HCl (pH8.0).

**Storage condition :** Stable at room temperature.  
For long time storage, – 20°C is recommended.

**Recommended usage:** Add 1/6 volume of 6× Loading Dye to the DNA sample, then apply to a well on agarose gel.

DNA sample	5 $\mu$ L
6× BPB loading Dye	1 $\mu$ L
<hr/>	
Mix	
↓	
Load to a well.	

Ver. 1.4