

Product Name:	Alkaline Phosphatase	(PAP)
	(from Shewanella sp	v. SIB1)*

Code No:	DE110		
Size:	1,000 unit		
<b>Concentration:</b>	5unit/µl		
Supplied reagent	<ul> <li>10×Alkaline Phosphatase Reaction Buffer</li> </ul>		
	Dilution Buffer		
Storage:	Store at -20°C		
Source:	<i>E.coli</i> harboring the plasmid encoding the gene of alkaline phosphatase from a psychrotrophic bacterium <i>Shewanella sp.</i> SIB1 (PAP).		
Unit Definition:	One unit is defined as the amount required to hydrolyzed 1.0 µmole p-nitrophenyl phosphate per 1 minute in glycine-NaOH buffer at pH10.5 and 37°C.		
Assay conditions:	The reaction mixture (100µl) contains 50mM glycine-NaOH buffer, pH 10.5, 5mM MgCl <sub>2</sub> , 0.5mM ZnCl <sub>2</sub> , 100mM KCl, 5mM p-nitrophanyl phosphate.		
Storage Buffer:	10mM Tris-HCl pH7.5 0.025mM ZnCl <sub>2</sub>		

#### **Contaminants:**

**Dnase:** When 0.5µg of  $\lambda$ /*Hind* III digest was incubated with 10units of this enzyme in a 40µl reaction mixture for 18 hours at 37°C, no degradation of the DNA fragment is observed on agarose gel electrophoresis.

**Rnase:** No RNase activity is observed by the use of RnaseAlert assay (Ambion). In this assay the reaction mixture containing the fluorescent-labeled RNA substrate was incubated with 10units of this enzyme for 1 hours at  $37^{\circ}$ C.

# \* Licensed Under Japan Patent NO. 2001-172653

0.25mM

50%

MgCl<sub>2</sub> glycerol



### **Composition of Supplied Reagent:**

 $10 \times Alkaline Phosphatase Reaction Buffer (Store at -20°C)$ 

 1. 5M
 Tris-HCl, pH7.3

 125mM
 glycine

 0. 5%
 TritonX-100

 0.25mM
 ZnCl2

 2.5mM
 MgCl2

 60mM
 NiCl2

# Dilution Buffer (1×Reaction Buffer, Store at -20°C)

12.5mM     glycine       0.05%     TritonX-100       0.025mM     ZnCl2       0.25mM     MgCl2       6mM     NiCl2	150mM	Tris-HCl, pH7.3
0.05%         TritonX-100           0.025mM         ZnCl2           0.25mM         MgCl2           6mM         NiCl2	12.5mM	glycine
0.025mM         ZnCl2           0.25mM         MgCl2           6mM         NiCl2	0.05%	TritonX-100
0.25mM MgCl <sub>2</sub> 6mM NiCl <sub>2</sub>	0.025mM	$ZnCl_2$
6mM NiCl <sub>2</sub>	0.25mM	MgCl <sub>2</sub>
	6mM	NiCl <sub>2</sub>



### Kit Manual

# **Dephosphorylation of 5' end by PAP:**

Linearized plasmid DNA \*1 88µl (5' termini 10pmol)\*2

10×Reaction buffer 10μl

Diluted PAP solution<sup>\*3</sup>  $2\mu l^{*4}$ 

Incubate for 30 min at 60°C (blunt end and 3'-protruding end)

or at 37°C (5'-protruding end)

Inactivation by incubation for 5min at  $95^{\circ}C^{*5}$ 

Dephosphorylated Plasmid DNA

# **Ligation Reaction and Transformation:**

	Dephosph	orylated Plasmid I	DNA	1μl (0.05pmol)* <sup>6</sup>
	Insert DNA (foreign)			1μl (0.05~0.3pmol)* <sup>7</sup>
	10×Ligation buffer* <sup>8</sup>			1µl
	T4 DNA li	gase		0.5 Weiss unit
	10mM AT	Р		1µl
	$H_2O$		up to	10µl
		Incubate for 1 ho	ure at 12°C	(blunt end)
		or at 16°C (5'-pro	otruding en	d and 3'-protruding end)
	2µl (0.0	lpmol) <mark>*</mark> <sup>9</sup>		
		Add 100µl of cor	npetent cel	1
		Store on ice for 3	0min	
		Heat shock for 30	Osec at 42°C	C
		Store on ice 1 min	ı	
		Add 0.9ml of SO	C medium	
		Incubate for 1hr	at 37°C	
		Spread on a selec	tive plate	
It	ncubate ove	rnight at 37°C		

\*1 Before the dephosphorylation, the complete digestion of the plasmid DNA should be confirmed by the agarose gel electrophoresis. Actually, it is difficult to completely digest vector DNA less than 4 kb in size. We recommend to purify the linearized vector DNA from agarose gel after electrophoresis by using Gel



example,

should

be

#### Indicator (Code No. DM510).

Restriction enzyme buffer (such as low buffer, medium buffer and high buffer) and 1×TE buffer is permissible as a buffer solution of the linearized plasmid DNA.

\*2 Amount and length of the linearized plasmid DNA

Table 1			
1kb linear DNA	10pmol of 5'-tremini = $3.3\mu g$		
2kb linear DNA	10pmol of 5'-tremini = $6.6\mu g$		
3kb linear DNA	10pmol of 5'-tremini = $9.9\mu g$		
4kb linear DNA	10pmol of 5'-tremini = $13.2 \mu g$		

For

10pmole of 5'-termini of the linearized pUC18 (2.69kb) is 8.8µg.

### \*3 Amount of PAP

Amount of PAP should be modified depending on the kind of termini and the amount of 5'-termini of linearized plasmid DNA. The following amounts are recommended:

Table 2				
	Units of PAP			
(10pmol)	1.0 units (37°C 30min)			
(10pmol)	2.5 units (60°C 30min)			
(10pmol)	5. 0 units (60°C 30min)			
	(10pmol) (10pmol) (10pmol)			

# PAP

diluted with dilution buffer (1×reaction buffer) according to the table above. As the amount of PAP described in the table above are about ten times as much as that of the minimum effective amount, the condition is sufficient for complete dephosphorylation. If you use over ten times amount of PAP shown in the table 2, we recommend phenol extraction to inactivate PAP completely instead of heat inactivation (see \*4 and \*5).

\*4 Non-diluted PAP solution should be added not to exceed 10% in a volume of the final reaction buffer. Glycerol in non-diluted PAP solution and high concentration of PAP protein hamper heat inactivation of PAP. If the reaction mixture contains non-diluted PAP solution up to 20% of its volume, the activity of about  $1/7,000 \sim 1/50,000$  still remains. When you use non-diluted PAP solution over 10% volume of reaction mixture, we recommend phenol extraction for complete inactivation of PAP (see \*5).

# \*5 Heat inactivation

At least 10  $\mu$ l of reaction mixture is required in 0.5ml tube for heat inactivation. Before heat inactivation, you should spin down solution in the tube. After heat inactivation, the white precipitate of inactivated PAP protein is observed. It is not necessary to be removed, because the precipitate does not affect the next ligation reaction.



Nearly all of PAP are inactivated by incubation at  $95^{\circ}$ C for 5min after dephosphorylation and trace amount activity,  $1/100,000 \sim 1/400,000$ , remains. It is negligible level in the next ligation reaction, because the residual PAP activity can only dephosphorylate insert DNA (foreign) less than 1/10,000 of the initial vector plasmid DNA in the molar amount, if you follow the above protocol.

If you use over ten times amount of enzyme the condition described in table 2 or use over 10% volume of final reaction mixture, we recommend phenol extraction to inactivate PAP completely instead of heat inactivation as follows;

Phenol Extraction:
Alkaline phosphatase reacton product
Phenol/chloroform (1:1) extraction (2 times)
Chloroform extraction (3 times)
Add 10ul 3M Na Acetate pH5.2 or 10ul 1M NaCl
Add 70µl isopropanol
Centrifuge at 14,500 rpm for 10min
Wash with 70% ethanol
Dry under vacuum
Dephosphorylated Plasmid DNA

PAP is easily and completely inactivated even by one time phenol extraction. On the other hand, BAP is known to be much more resistant to phenol extraction as well as heat inactivation. For example, about 2.5% of BAP remain active after one time phenol extraction and following ethanol precipitation (http://www.biodynamics.co.jp/).

\*6 It is possible that alkaline phosphatase reaction mixture is directly added to ligation reaction mixture up to 30% for cohesive end ligation, up to 10% for blunt end ligation of the final reaction volume without interference.

The commercially available ligation kit such as <sup>DynaExpress</sup>DNA Ligation Kit ver.2 (BioDynamics Laboratory.Inc, DS110), DNA Ligation Kit ver.2 (TaKaRa), Ligation High (TOYOBO CO., LTD), Quick Ligation<sup>TM</sup> Kit (New England Biolabs. Inc) and LigaFast<sup>TM</sup> Rapid DNA Ligation Kit (Promega) can be effectively used in instead of this ligation protocol. During you use these kits, the orange precipitate derived from reducted Ni is observed. It is not necessary to remove the precipitate because it does not affect the ligation and the next transformation process.

\*7 Chill the reaction mixture of insert DNA and vector DNA to  $0^{\circ}$ C before ligation procedure.

The molar ratio of vector DNA to insert DNA should be between 1:1 and 1:6. The final DNA concentration of vector DNA and insert DNA should be between  $1ng/\mu l$  and  $10ng/\mu l$  for an effective ligation.

LigaFast<sup>TM</sup> and Quick Ligation<sup>TM</sup> are trademarks of Promega Corporation and New England Biolabs, Inc., respectively.

\*8 10×Ligation buffer

660mM Tris-HCl (pH7.6), 66mM MgCl<sub>2</sub>, 100mM DTT



\*9 Amount of ligation mixture should be added in a volume not to exceed 10% of that of competent cells. If you carry out transformation of *E.coli* by electroporation, we recommend spin column purification of ligation mixture.