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Product Name: Japanese Cedar Pollen Allergen Cry j 2, Purified

Code No.: HBL-C-2

Size: 25 µg, 250 µl/vial (Frozen Liquid)

This product is for research use only

Specification and Use

Preparation & form :

Cry j 2 antigen was extracted from Japanese cedar (*Cryptomeria japonica*) pollen and purified by ammonium sulfate precipitation, ion-exchange chromatography and anti-Cry j 2 monoclonal antibody immobilized column chromatography.

Purified Cry j 2 allergen was prepared in the concentration of 100 µg/ml with phosphate-buffered saline and frozen. No stabilizers or preservatives are contained.

Molecular weight :

37 kDa, as determined by SDS-PAGE under a non-reduced condition.

Purity:

Greater than 95%.

Cry j 1 antigen is contained less than 0.1%, measured with ELISA.

Storage :

Store below -20°C.

Avoid multiple freeze/thaw cycles by storage in appropriate aliquots.

This material is unstable at 5°C.

References :

- Sakaguchi, M et al.: Fluorometric enzyme-linked immunosorbent assay for the measurement of IgE antibody to *Cryptomeria japonica* pollen. [Article in Japanese] *Arerugi*, 34(4):233-237, 1986
- Sawatani, M et al. : Immunological and physiological properties of Cry j II, the second major allergen of Japanese cedar pollen. [Article in Japanese] *Arerugi*, 42(6):738-747, 1993
- Sone, T et al.: Cloning and sequencing of cDNA coding for Cry j I, a major allergen of Japanese cedar pollen. *Biochem, Biophys. Res. Comm.* 199(2):619-625, 1995
- Namba, M. et al.: Molecular cloning of the second major allergen, Cry j II, from Japanese cedar pollen. *FEBS Letter*, 353(2):124-128, 1994

Caution ! :

This material is research use only, not for clinical or diagnostic use.

This material has allergenicity.