

PRODUCT INFORMATION

Product Name: RNase-free Water, non DEPC

This product is research use only

Code No and Size:

Code	Product Name	Size
DR120	RNase-free Water, non DEPC, 1 ml × 5	1 ml tube × 5, total 5 ml
DR125	RNase-free Water, non DEPC, 50 ml × 2	50 ml bottle × 2, total 100 ml
DR127	RNase-free Water, non DEPC, 500 ml	500 ml bottle × 1

Storage: Store at room temperature

Description:

RNase-free Water is suitable to experiment and research of RNA. The RNase-free Water, non DEPC, was prepared from highly purified water by ultrafiltration without DEPC treatment. The product was confirmed to be free from RNase and DNase. It will be stable for at least one year from the date of receipt if stored in unopened packing.

Quality Control:

- DNase activity: None DNase activity was detected after incubation plasmid DNA and this product for 18 hr at room temperature.
- RNase activity: None RNase activity was detected after incubation ribosomal RNA and this product for 18 hr at room temperature.

Ver. 1.2